

Product Datasheet

Product:	Alpha 1 Acid Glycoprotein, Cynomolgus Monkey
Catalog #:	16-16-010700-CM
Form:	Lyophilized
Available Packaging:	100 µg. Larger aliquots are available upon request.
Buffer:	Lyophilized from deionized water
Purity:	≥ 95% by SDS-PAGE
Molecular Weight:	55,000 Da
Protein Determination:	Extinction Coefficient (E) = 0.89 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Monkey Plasma
Testing:	Monkeys were individually tested for B-Virus (Herpes Virus Simiae) and found to be negative.

Product Link: <https://www.athensresearch.com/products/all/alpha-1-acid-glycoprotein-cynomolgus-monkey-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg