


Gel Scan of Representative Lot Transferrin, Rat

ART Catalog No. 16-16-032001-rat



SDS-PAGE

4-12% Bis-Tris gel, 1x MES

1. Transferrin, Rat 5 µg (Heated/Reduced)
2. Transferrin, Rat 10 µg (Heated/Reduced)
3. Transferrin, Rat 20 µg (Heated/Reduced)
4. Standard
5. Blank
6. Transferrin, Rat 5 µg (Not Heated/Non-Reduced)
7. Transferrin, Rat 10 µg (Not Heated/Non-Reduced)
8. Transferrin, Rat 20 µg (Not Heated/Non-Reduced)

Protein Determination:

Biuret Method

Considerations:

Storage Conditions: ≤ -20 °C

Stability: > 1 year

Molecular Weight:

80,000 Da

Buffer:

20 mM sodium phosphate, pH 7.4, 150 mM NaCl

Physical Specifications:

Form: Lyophilized

Purity: ≥ 98% by SDS-PAGE.

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg

Product Datasheet

Product:	Transferrin, Rat
Catalog #:	16-16-032001-rat
Form:	Lyophilized
Available Packaging:	1 mg. Large aliquots are available upon request
Buffer:	20 mM Na ₂ HPO ₄ , pH 7.4, 150 mM NaCl
Purity:	≥ 98% by SDS-PAGE.
Molecular Weight:	80,000 Da
Protein Determination:	Biuret Method
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Rat Plasma
Testing:	Inquire for available testing documentation for this material.

Product Link: <https://www.athensresearch.com/products/all/transferrin-rat-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg