

Product Datasheet

Product:	Hemopexin, Rat
Catalog #:	16-16-080513-Rat
Form:	Lyophilized
Available Packaging:	100 and 500 µg and 1 mg. Larger aliquots are available upon request.
Buffer:	Lyophilized from 11.9 mM sodium phosphate, 137 mM NaCl, 2.7 mM KCl, pH 7.4.
Purity:	≥ 95% by SDS-PAGE.
Molecular Weight:	57,000 Da
Protein Determination:	Extinction Coefficient (E) = 1.97 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Rat Serum
Testing:	Inquire for available testing documentation for this material.

Product Link: <https://www.athensresearch.com/products/all/hemopexin-rat-serum>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg