


Gel Scan of Representative Lot Immunoglobulin G3 (IgG3)

ART Catalog No. 16-16-090707-3



SDS-PAGE

4-12% Bis-Tris gel, 1x MES

1. Immunoglobulin G3 (IgG3) 5 µg (Heated/Reduced)
2. Immunoglobulin G3 (IgG3) 10 µg (Heated/Reduced)
3. Immunoglobulin G3 (IgG3) 20 µg (Heated/Reduced)
4. Standard
5. Blank
6. Immunoglobulin G3 (IgG3) 5 µg (Not Heated/Non-Reduced)
7. Immunoglobulin G3 (IgG3) 10 µg (Not Heated/Non-Reduced)
8. Immunoglobulin G3 (IgG3) 20 µg (Not Heated/Non-Reduced)

Protein Determination:

Extinction Coefficient (E) = 1.36 (0.1% at 280 nm, 1 cm pathway)

Molecular Weight:

170,000 Da

Physical Specifications:

Form: Liquid

Purity: ≥ 95% by SDS-PAGE

Considerations:

Storage Conditions: 2 - 8 °C, Do not freeze

Stability: > 1 year

Buffer:

Liquid in 20 mM sodium phosphate, pH 7.4, with 150 mM NaCl and 0.05% sodium azide.

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg

Product Datasheet

Product:	Immunoglobulin G3 (IgG3)
Catalog #:	16-16-090707-3
Form:	Liquid
Available Packaging:	100 µg and 1 mg. Larger aliquots are available upon request.
Buffer:	Liquid in 20 mM sodium phosphate, pH 7.4, with 150 mM NaCl and 0.05% NaN ₃ .
Purity:	≥ 95% by SDS-PAGE
Molecular Weight:	170,000 Da
Protein Determination:	Extinction Coefficient (E) = 1.36 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	2 - 8 °C, Do not freeze
Stability:	> 1 year
Source:	Human Plasma
Testing:	Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request.

Product Link: <https://www.athensresearch.com/products/all/immunoglobulin-g3-human-plasma-igg3>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg