


Gel Scan of Representative Lot Immunoglobulin G (IgG), Fab '2

ART Catalog No. 16-16-090707-FAB2



SDS-PAGE

4-12% Bis-Tris gel, 1x MES

1. Immunoglobulin G (IgG), Fab '2 2 µg (Heated/Reduced)
2. Immunoglobulin G (IgG), Fab '2 5 µg (Heated/Reduced)
3. Immunoglobulin G (IgG), Fab '2 10 µg (Heated/Reduced)
4. Immunoglobulin G (IgG), Fab '2 20 µg (Heated/Reduced)
5. Standard
6. Blank
7. Immunoglobulin G (IgG), Fab '2 2 µg (Not Heated/Non-Reduced)
8. Immunoglobulin G (IgG), Fab '2 5 µg (Not Heated/Non-Reduced)
9. Immunoglobulin G (IgG), Fab '2 10 µg (Not Heated/Non-Reduced)
10. Immunoglobulin G (IgG), Fab '2 20 µg (Not Heated/Non-Reduced)

Protein Determination:

Extinction Coefficient (E) = 1.36 (0.1% at 280 nm,
1 cm pathway)

Considerations:

Storage Conditions: ≤ -20 °C
Stability: > 5 years

Molecular Weight:

110,000 Da

Buffer:

Frozen in 20 mM sodium phosphate, pH 7.4, with
150 mM NaCl and 0.05% sodium azide.

Physical Specifications:

Form: Frozen

Purity: ≥ 95% by SDS-PAGE. Single arc by IEP
against antisera to whole human serum.

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg

Product Datasheet

Product:	Immunoglobulin G (IgG), Fab '2
Catalog #:	16-16-090707-FAB2
Form:	Frozen
Available Packaging:	1 mg. Larger aliquots are available upon request.
Buffer:	Frozen in 20 mM sodium phosphate, pH 7.4, with 150 mM NaCl and 0.05% NaN ₃ .
Purity:	≥ 95% by SDS-PAGE. Single arc by IEP against antisera to whole human serum.
Molecular Weight:	110,000 Da
Protein Determination:	Extinction Coefficient (E) = 1.36 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	≤ -20 °C
Stability:	> 5 years
Source:	Human Plasma
Testing:	Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request.

Product Link: <https://www.athensresearch.com/products/all/immunoglobulin-g-fab-2-human-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg