


Gel Scan of Representative Lot

Immunoglobulin G (IgG), FC Fragment, Immunogen Grade

ART Catalog No. 16-16-090707-FC-IMM



SDS-PAGE

4-12% Bis-Tris gel

1. Immunoglobulin G (IgG), FC Fragment, Immunogen Grade 5 µg (Heated/Reduced)
2. Immunoglobulin G (IgG), FC Fragment, Immunogen Grade 10 µg (Heated/Reduced)
3. Immunoglobulin G (IgG), FC Fragment, Immunogen Grade 20 µg (Heated/Reduced)
4. Standard
5. Blank
6. Immunoglobulin G (IgG), FC Fragment, Immunogen Grade 5 µg (Not Heated/Non-Reduced)
7. Immunoglobulin G (IgG), FC Fragment, Immunogen Grade 10 µg (Not Heated/Non-Reduced)
8. Immunoglobulin G (IgG), FC Fragment, Immunogen Grade 20 µg (Not Heated/Reduced)

Protein Determination:

Extinction Coefficient (E) = 1.35 (0.1% at 280 nm, 1 cm pathway)

Considerations:

Storage Conditions: ≤ -20 °C

Stability: > 1 year

Molecular Weight:

50,000 Da

Buffer:

Frozen in 50 mM Tris-HCl, pH 8.0, with 200 mM NaCl and 0.05% sodium azide.

Physical Specifications:

Form: Frozen

Purity: 99% by SDS-PAGE; Single arc by IEP;

Immunogen Grade per Western Blot

Your Distributor

ORDERING
 Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
 PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
 Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
 Oberaltenallee 8
 D-22081 Hamburg
 HRB 98170 Amtsgericht Hamburg

Product Datasheet

Product:	Immunoglobulin G (IgG), FC Fragment, Immunogen Grade
Catalog #:	16-16-090707-FC-IMM
Form:	Frozen
Available Packaging:	1, 5, and 10 mg. Larger aliquots are available upon request.
Buffer:	Frozen in 50 mM Tris-HCl, pH 8.0, with 200 mM NaCl and 0.05% NaN ₃ .
Purity:	99% by SDS-PAGE; Single arc by IEP; Immunogen Grade per Western Blot
Molecular Weight:	50,000 Da
Protein Determination:	Extinction Coefficient (E) = 1.35 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Human Plasma
Testing:	Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request.

Product Link: <https://www.athensresearch.com/products/all/immunoglobulin-g-fc-fragment-human-plasma-immunogen-grade>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg