

Product Datasheet

Product:	Immunoglobulin G (IgG), Mouse
Catalog #:	16-16-090707-mouse
Form:	Lyophilized
Available Packaging:	1 and 5 mg. Larger aliquots are available upon request.
Buffer:	20 mM Na ₂ HPO ₄ , pH 7.4 with 150 mM NaCl and 0.05% NaN ₃
Purity:	≥ 95% by SDS-PAGE
Molecular Weight:	150,000 Da
Protein Determination:	Extinction Coefficient (E) = 1.36 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Mouse Plasma
Testing:	Inquire for available testing documentation for this material.

Product Link: <https://www.athensresearch.com/products/all/immunoglobulin-g-mouse-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg