


Gel Scan of Representative Lot Kininogen, LMW

ART Catalog No. 16-16-110914-L



SDS-PAGE

4-12% Bis-Tris gel, 1x MES

1. Kininogen, LMW 5 µg (Heated/Reduced)
2. Kininogen, LMW 10 µg (Heated/Reduced)
3. Kininogen, LMW 20 µg (Heated/Reduced)
4. Standard
5. Blank
6. Kininogen, LMW 5 µg (Not Heated/Non-Reduced)
7. Kininogen, LMW 10 µg (Not Heated/Non-Reduced)
8. Kininogen, LMW 20 µg (Not Heated/Non-Reduced)

Protein Determination:

Extinction Coefficient (E) = 0.59 (0.1% at 280 nm,
1 cm pathway)

Molecular Weight:

65,000 Da

Physical Specifications:

Form: Lyophilized

Purity: ≥ 95% by SDS-PAGE

Considerations:

Storage Conditions: ≤ -20 °C

Stability: > 1 year

Buffer:

Lyophilized from 10 mM sodium acetate, pH 5.5,
with 200 mM NaCl.

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg

Product Datasheet

Product:	Kininogen, LMW
Catalog #:	16-16-110914-L
Form:	Lyophilized
Available Packaging:	100 µg and 1 mg. Larger aliquots are available upon request.
Buffer:	Lyophilized from 10 mM sodium acetate, pH 5.5, with 200 mM NaCl.
Purity:	≥ 95% by SDS-PAGE
Molecular Weight:	65,000 Da
Protein Determination:	Extinction Coefficient (E) = 0.59 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Human Plasma
Testing:	Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request.

Product Link: <https://www.athensresearch.com/products/all/kininogen-lmw-human-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg