

Product Datasheet

Product:	Apolipoprotein CI
Catalog #:	16-16-120301
Form:	Lyophilized
Available Packaging:	100 µg and 1 mg. Larger aliquots are available upon request.
Buffer:	Lyophilized from 10 mM NH ₄ HCO ₃ , pH 7.4
Purity:	≥ 95% by SDS-PAGE
Molecular Weight:	6,600 Da
Protein Determination:	Modified Lowry based on BSA Standard
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Human Plasma
Testing:	Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request.

Product Link: <https://www.athensresearch.com/products/all/apolipoprotein-ci-human-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg