

Product Datasheet

Product:	Apolipoprotein E
Catalog #:	16-16-120500
Form:	Frozen
Available Packaging:	50 and 500 µg and 1 mg. Larger aliquots are available upon request.
Buffer:	Frozen in 50 mM NH ₄ HCO ₃ , pH 8.0
Purity:	≥ 95% by SDS-PAGE
Molecular Weight:	34,000 Da
Protein Determination:	Extinction Coefficient (E) = 1.31 (0.1% at 280 nm, 1 cm pathway)
Storage/Handling:	≤ -70 °C
Stability:	> 1 year
Source:	Human Plasma
Testing:	Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request.

Product Link: <https://www.athensresearch.com/products/all/apolipoprotein-e-human-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg