

Product Datasheet

| | |
|-------------------------------|---|
| Product: | Rheumatoid Factor IgM |
| Catalog #: | 16-16-180603 |
| Form: | Liquid |
| Available Packaging: | 250 µL and 1 mL. Larger aliquots are available upon request. |
| Buffer: | Liquid in 20 mM sodium phosphate, pH 7.4, with 150 mM NaCl and 0.09% NaN ₃ |
| Molecular Weight: | 970,000 Da |
| Protein Determination: | Extinction Coefficient (E) = 1.18 (0.1% at 280 nm, 1 cm pathway) |
| Storage/Handling: | 2 - 6 °C |
| Stability: | > 1 year |
| Source: | Human Plasma |
| Testing: | Donor material is obtained from suppliers that perform a comprehensive infectious disease screening panel. Each unit of plasma is tested by the supplier and found non-reactive for HIV-1/2, Hepatitis B surface antigen (HBsAg), Hepatitis C (HCV), Hepatitis A (HAV), and syphilis (RPR). Documentation of testing is maintained and is available upon request. |

Product Link: <https://www.athensresearch.com/products/all/rheumatoid-factor-igm-human-plasma>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg