

Product Datasheet

Product:	S-100b Protein, Bovine
Catalog #:	16-BR-190802-bovine
Form:	Lyophilized
Available Packaging:	100 µg and 1 mg. Larger aliquots are available on request.
Buffer:	50 mM Tris-HCl, pH 7.5; 50 mM NaCl; 100 µM DTT
Purity:	≥ 95% by SDS-PAGE
Molecular Weight:	6,000 Da
Protein Determination:	Bradford Assay
Storage/Handling:	≤ -20 °C
Stability:	> 1 year
Source:	Bovine Brain
Testing:	Prepared from tissue collected in the contiguous 48 United States from animals under 30 months of age at a USDA-inspected facility. Animals were inspected ante mortem and found free from obvious signs of infectious and contagious diseases.

Product Link: <https://www.athensresearch.com/products/all/s-100b-protein-bovine-brain>

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Your Distributor

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de



TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

ACTIVE BIOSCIENCE
PROTEINS · CELLS · ANTIBODIES

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg