



Human JNK2/SAPK1

Synonyms: Mitogen-activated protein kinase 9, EC 2.7.11.24, Stress-activated protein kinase JNK2, c-Jun N-terminal kinase 2, JNK-55, JNK2, JNK2A, JNK2B, PRKM9, JNK2BETA, p54aSAPK, JNK2ALPHA, MAPK9, SAPK1.

PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING

Size	Order#	Lot#	Expiry Date
1 µg	4029.950.001		
5 µg	4029.950.005		
200 µg	4029.950.200		

Please enquire for bulk quantities and other vial sizes.

Description

SAPK1/JNK2 Human Recombinant (stress-activated protein kinase/c-jun kinase) is a non-glycosylated polypeptide having a molecular mass of 49.5 kDa.

- **Source** *E. coli*
- **Purity** ≥ 95 % (SDS-PAGE)
- **Buffer** JNK2/SAPK1 (1.3 mg/ml) is supplied in 50 mM Tris-HCl, 100 mM NaCl, 1 mM DTT, 50 % glycerol, pH 8*
- **Physical State** Sterile filtered clear solution

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. **Please avoid repeated freeze-thaw cycles.**

Usage: For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.