



## Anti-Human KDR / VEGFR-2, bt (#3) (MAB)

Synonyms: vascular endothelial growth factor receptor-2 , KDR, FLK1, CD309, VEGF receptor 2, VEGFR2, kinase insert domain protein receptor

*PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING*

| Size                                                     | Order#       | Lot# | Expiry Date |
|----------------------------------------------------------|--------------|------|-------------|
| 50 µg                                                    | 1620.551.050 |      |             |
| Please enquire for bulk quantities and other vial sizes. |              |      |             |

### Description

VEGF R1 (Flt-1), VEGF R2 (KDR/Flk-1), and VEGF R3 (Flt-4) belong to the class III subfamily of receptor tyrosine kinases (RTKs). All three receptors contain seven immunoglobulin-like repeats in their extracellular domain and kinase insert domains in their intracellular region. They are best known for regulating VEGF family-mediated vasculogenesis, angiogenesis, and lymphangiogenesis. They are also mediators of neurotrophic activity and regulators of hematopoietic development. Human VEGF R2 is thought to be the primary inducer of VEGF-mediated blood vessel growth, while VEGF R3 plays a significant role in VEGF-C and VEGF-D-mediated lymphangiogenesis.

- **Source** Mouse
- **Isotype** IgG1
- **Clone** AB-45 I4

### Biological Activity

FACS: Use at 2-10 µg/ml; Western Blot: Use at 2-5 µg/ml

### Reconstitution

Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.1-1.0 mg/ml.

**Usage:** For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

---

\*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING  
Tel.: +49 40 4320 8448 - 0  
order@active-bioscience.de  
www.active-bioscience.de

TECHNICAL SUPPORT  
Tel.: +49 40 4320 8448 - 11  
support@active-bioscience.de

Active Bioscience GmbH  
Oberaltenallee 8  
D-22081 Hamburg  
HRB 98170 Amtsgericht Hamburg