



Anti-Human KDR / VEGFR-2 / scFV (single chain)

Synonyms: vascular endothelial growth factor receptor-2 , KDR, FLK1, CD309, VEGF receptor 2, kinase insert domain protein receptor

PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING

Size	Order#	Lot#	Expiry Date
100 µg	1619.856.100		
Please enquire for bulk quantities and other vial sizes.			

Description

Anti-human VEGFR-2 single chain (VEGFR-2 sc-Fv) antibody was constructed from cDNA of mice immunized with highly purified soluble VEGFR-2 protein expressed in insect cells. Anti-human VEGFR-2 sc-Fv was produced in *E. coli* and affinity purified.

- **Source** *E. coli*
- **Clone** AB-Sfdpncjobou!JH

Biological Activity

FACS: Use at 0.5-5µg/ml; ELISA

Usage: For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING
Tel.: +49 40 4320 8448 - 0
order@active-bioscience.de
www.active-bioscience.de

TECHNICAL SUPPORT
Tel.: +49 40 4320 8448 - 11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg