



## Anti-Murine CD105 / Endoglin, soluble PE

Synonyms: Eng endoglin, CD105, AI528660, AI662476, S-endoglin, FITC

*PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING*

Size	Order#	Lot#	Expiry Date
25 µg	1260.662.025		
50 µg	1260.662.050		

Please enquire for bulk quantities and other vial sizes.

### Description

The MJ7/18 monoclonal antibody reacts with the mouse CD105 molecule, also known as Endoglin. This 90 kDa disulfide-linked homodimer is expressed by vascular endothelial cells. It is suggested that CD105 functions in adhesion and embryonic angiogenesis.

- **Source** Rat
- **Isotype** IgG2 a
- **Clone** AB-NK8029

**Usage:** For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

---

\*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING  
Tel.: +49 40 4320 8448 - 0  
order@active-bioscience.de  
www.active-bioscience.de

TECHNICAL SUPPORT  
Tel.: +49 40 4320 8448 - 11  
support@active-bioscience.de

Active Bioscience GmbH  
Oberaltenallee 8  
D-22081 Hamburg  
HRB 98170 Amtsgericht Hamburg