

Equine Leptin

Synonyms: OB Protein, Obesity Protein, OBS, Obesity factor.

PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING

Size	Order #	Lot #	Expiry Date
20 µg	1650.900.020		
100 µg	1650.900.100		
1 mg	1650.900.199		

Please enquire for bulk quantities and other vial sizes

Description

Leptin Horse Recombinant produced in *E. Coli* is a single, non-glycosylated, polypeptide chain containing 146 amino acids and having a molecular mass of 16 kDa.

- **Source** *E. Coli*
- **Purity** ≥ 95 % (SDS-PAGE, SEC-HPLC)
- **Endotoxin level** ≤ 0.1ng/µg (≤ 1EU/µg)
- **Buffer** The protein was lyophilized from a concentrated (1mg/ml) solution with 0.0045mM NaHCO₃*
- **Physical state** Sterile filtered, lyophilized

Biological Activity

Biological active as evidenced by inducing proliferation of BAF/3 cells stably transfected with the long form of human leptin receptor.

Reconstitution

We recommend a quick spin followed by reconstitution in sterile water to a concentration of at least 100 µg/ml, which can then be further diluted to other aqueous solutions. Do not vortex.

Stability

Lyophilized Leptin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Leptin should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). **Please avoid repeated freeze-thaw cycles.**

Amino Acid Sequence

The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Val-Pro-Ile-Gln.

Usage: For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING
Tel.: +49 40 43208448-0
order@active-bioscience.de
www.active-bioscience.de

TECHNICAL SUPPORT
Tel.: +49 40 43208448-11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg