

Human Growth Hormone Placental 22 kDa

Synonyms: GHL, GHV, GH-V, hGH-V, PGH.

PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING

Size	Order #	Lot #	Expiry Date
10 µg	1448.950.010		
50 µg	1448.950.050		
1 mg	1448.950.199		

Please enquire for bulk quantities and other vial sizes

Description

Placental HGH 22kDa Recombinant produced in *E. Coli* is a single, non-glycosylated, polypeptide chain containing 192 amino acids and having a molecular mass of 22367 Dalton. Predicted pI=7.80. Placental Growth Hormone has diminished lactogenic (prolactin receptor mediated) activity characteristic to pituitary GHs.

- **Source** *E. Coli*
- **Purity** ≥ 98 % (SDS-PAGE, SEC-HPLC)
- **Endotoxin level** ≤ 0.1ng/µg (≤ 1EU/µg)
- **Buffer** The protein was lyophilized from a concentrated (1mg/ml) solution with 0.0045mM NaHCO₃ previously adjusted pH 8-9*
- **Physical state** Sterile filtered, lyophilized

Reconstitution

It is recommended to reconstitute the lyophilized Placental HGH in 0.4% NaHCO₃ or water adjusted to pH 9, not less than 100 µg/ml, which can then be further diluted to other aqueous solutions.

Stability

Lyophilized GH placental although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution and filter sterilization HGH placental can be stored at 4°C for up to 4 weeks. For long term storage and more diluted solutions it is recommended to add a carrier protein (0.1% HSA or BSA). **Please avoid repeated freeze-thaw cycles.**

Amino Acid Sequence

The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Phe-Pro-Thr-Ile.

Usage: For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING
Tel.: +49 40 43208448-0
order@active-bioscience.de
www.active-bioscience.de

TECHNICAL SUPPORT
Tel.: +49 40 43208448-11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg