

## Human SCF, cct-premium

Synonyms: KITLG, SF, MGF, SCF, FPH2, KL-1, Kitl, SHEP7, kit-ligand, stem cell factor, Mast cell growth factor, c-Kit ligand

*PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING*

Size	Order #	Lot #	Expiry Date
2 µg	1884.950.002		
10 µg	1884.950.010		
50 µg	1884.950.050		
100 µg	1884.950.100		
500 µg	1884.950.500		
1 mg	1884.950.199		

Please enquire for bulk quantities and other vial sizes

### Description

Stem cell factor (SCF), also known as c-kit ligand (KL), mast cell growth factor (MGF), and steel factor (SLF), is a widely expressed 28-40 kDa type I transmembrane glycoprotein. It promotes the survival, differentiation, and mobilization of multiple cell types including myeloid, erythroid, megakaryocytic, lymphoid, germ cell, and melanocyte progenitors. SCF is a primary growth and activation factor for mast cells and eosinophils. Mature mouse SCF consists of a 189 amino acids (aa) extracellular domain (ECD), a 23 aa transmembrane segment, and a 36 aa cytoplasmic tail. The ECD shows both N-linked and O-linked glycosylation. Proteolytic cleavage at two alternate sites in the extracellular juxtamembrane region releases a 25 kDa soluble molecule which is comparable to the only form produced by Steel-dickie mutant mice. An alternatively spliced isoform of mouse SCF lacks 28 aa that encompasses the primary proteolytic recognition site. Within the ECD of the short isoform (corresponding to this recombinant protein), mouse SCF shares 93% aa sequence identity with rat SCF and 72% to 75% with canine, feline, and human SCF. Rat SCF is active on mouse and human cells, but human SCF is only weakly active on mouse cells. Non-covalent dimers of transmembrane or soluble SCF interact with the receptor tyrosine kinase SCF R/c-kit to trigger receptor dimerization and signaling. SCF assists in the recovery of cardiac function following myocardial infarction by increasing the number of cardiomyocytes and vascular channels.

- **Source** *E. Coli*
- **Purity** ≥ 98 % (SDS-PAGE, silver stained)
- **Endotoxin level** ≤ 0.01ng/µg (≤ 0.1EU/µg)

### Biological Activity

The ED<sub>50</sub> as determined by the dose-dependent stimulation of human TF-1 cells is < 2ng/ml, corresponding to a specific activity of 5x10<sup>5</sup> IU/mg.

### Reconstitution

Centrifuge vial prior to opening. Human SCF should be reconstituted in water to a concentration of 0.1 mg/ml. This solution can be further diluted in water or other buffer solutions or stored at -20 °C.

### Amino Acid Sequence

aa26t0aa15

**Usage:** For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

\*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING  
Tel.: +49 40 43208448-0  
order@active-bioscience.de  
www.active-bioscience.de

TECHNICAL SUPPORT  
Tel.: +49 40 43208448-11  
support@active-bioscience.de

Active Bioscience GmbH  
Oberaltenallee 8  
D-22081 Hamburg  
HRB 98170 Amtsgericht Hamburg