

Rat GRO-alpha / KC / CXCL1

Synonyms: Growth-regulated protein alpha, CXCL1, Melanoma growth stimulatory activity, MGSA, Neutrophil-activating protein 3, NAP-3, GRO-alpha(1-73), chemokine (C-X-C motif) ligand 1, GRO1, GROa, SCYB1, MGSA-a, MGSA alpha.

PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING

Size	Order #	Lot #	Expiry Date
5 µg	1455.970.005		
25 µg	1455.970.025		
1 mg	1455.970.199		

Please enquire for bulk quantities and other vial sizes

Description

CXCL1 Rat Recombinant produced in *E. Coli* is a single, non-glycosylated, polypeptide chain containing 73 amino acids and having a molecular mass of 7.8 kDa.

- **Source** *E. Coli*
- **Purity** ≥ 97 % (SDS-PAGE, RP-HPLC)
- **Endotoxin level** ≤ 0.1ng/µg (≤ 1EU/µg)
- **Buffer** Lyophilized from a 0.2µm filtered concentrated (1mg/ml) solution in 20mM PB, pH 7.4, and 150mM NaCl*
- **Physical state** Sterile filtered, lyophilized

Biological Activity

Determined by its ability to chemoattract rat neutrophils using a concentration range of 10-100ng/ml corresponding to a Specific Activity of 10,000-100,000IU/mg.

Reconstitution

We recommend a quick spin followed by reconstitution in sterile water to a concentration of at least 100 µg/ml, which can then be further diluted to other aqueous solutions. Do not vortex.

Stability

Lyophilized CXCL1 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution CXCL1 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). **Please avoid repeated freeze-thaw cycles.**

Amino Acid Sequence

APVANELRCQ CLQTVAGIHF KNIQSLKVMP PGPHTQTEV IATLKNAREA CLDPEAPMVQ KIVQKMLKGV
PK

Usage: For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING
Tel.: +49 40 43208448-0
order@active-bioscience.de
www.active-bioscience.de

TECHNICAL SUPPORT
Tel.: +49 40 43208448-11
support@active-bioscience.de

Active Bioscience GmbH
Oberaltenallee 8
D-22081 Hamburg
HRB 98170 Amtsgericht Hamburg