

## Rat Leptin Triple Antagonist

Synonyms:

*PLEASE NOTE: ALWAYS CENTRIFUGE VIAL BEFORE OPENING*

Size	Order #	Lot #	Expiry Date
10 µg	1649.970.010		
50 µg	1649.970.050		
1 mg	1649.970.199		

Please enquire for bulk quantities and other vial sizes

### Description

Leptin Antagonist Triple Mutant Rat Recombinant is a singly non-glycosylated polypeptide chain containing 146 amino and additional Ala at N-terminus acids and having a molecular mass of 16 kDa, Leptin was mutated, resulting in L39A/D40A/F41A mutant.

- **Source** *E. Coli*
- **Purity** ≥ 99 % (SDS-PAGE, Gel Filtration)
- **Endotoxin level** ≤ 0.1ng/µg (≤ 1EU/µg)
- **Buffer** The protein was lyophilized from a concentrated (0.85mg/ml) solution with 0.003mM NaHCO<sub>3</sub>\*
- **Physical state** White lyophilized (freeze-dried) powder

### Biological Activity

Leptin-induced proliferation of BAF/3 cells stably transfected with the long form of human Leptin receptor can be inhibited by recombinant rat Leptin Antagonist Triple Mutant. It also inhibits various Leptin effects in several in vitro bioassays.

### Reconstitution

It is recommended to reconstitute the lyophilized Leptin Antagonist Triple Mutant Rat Recombinant in sterile water or sterile 0.4% NaHCO<sub>3</sub> adjusted to pH 8-9, not less than 100µg/ml, which can then be further diluted with other aqueous solutions.

### Stability

Lyophilized Leptin-Antagonist Triple Mutant Rat Recombinant although stable at room temperature for several weeks, should be stored desiccated below -18°C. Upon reconstitution at > 0.1 Leptin mutant mg/ml and up to 2 mM and filter sterilization LEP-tA mutant can be stored at 4°C or even room temperature for several weeks making it suitable for long term infusion studies using osmotic pumps. At lower concentration addition of a carrier protein (0.1% HSA or BSA) is suggested. **Please avoid repeated freeze-thaw cycles.**

### Amino Acid Sequence

The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Val-Pro-Ile-Gln.

**Usage:** For research use only. Not for use in diagnostic or therapeutic procedures. Not for human use.

\*The Buffer may vary depending on the Lot #. Please contact our technical support if you have specific requirements.

ORDERING  
Tel.: +49 40 43208448-0  
order@active-bioscience.de  
www.active-bioscience.de

TECHNICAL SUPPORT  
Tel.: +49 40 43208448-11  
support@active-bioscience.de

Active Bioscience GmbH  
Oberaltenallee 8  
D-22081 Hamburg  
HRB 98170 Amtsgericht Hamburg